

# Quik-Tip™

Cone Sockets With XP Technology™

- **The Only Sharpen-Free Cone Socket Instruments On The Planet**

Some may claim to reduce the frequency of sharpening. XP Technology completely eliminates the need to resharpen.

- **Out-Performs Competitors Carbon Steel by up to 240 times!**

Results from an independent metallurgical testing lab found XP technology to be up to 240 times harder than our leading competitors carbon steel.

- **At A Fraction Of The Cost**

Replace only the tips that are worn-out at a fraction of the cost of an entirely new instrument.

- **Environmentally Friendly**

By replacing only the tips instead of the entire instrument, waste is greatly reduced.

- **Custom Combinations**

Create your perfect set-up with the most popular tip combinations with the Quik-Tip system.

- **Quik-Tip Cone Socket Installation**

Our Quik-Tip cone sockets offer easy installation that takes seconds, finger tighten and a 1/8th turn of the wrench and you're done.



See back for more information →

## The Only Cone Sockets Available With XP Technology

Creating The Only Sharpen-Free Cone Socket Instruments On The Planet.

- **Lasts Months Without Sharpening**

Have the same factory sharp edge every time you pick up your instrument.

- **More Comfort and Less Fatigue**

An XP instrument that remains super sharp can always be used with the same light pressure, when planing away calculus or tartar.

- **Same Special Care You Are Used To**


We recommend that our cone socket tips be cleaned and sterilized using the same methods you use on your XP instruments.

- **Save Up To \$4400 Per Year**

By eliminating the need for sharpening, some offices save up to \$4400 a year in labor costs.

- **Hartzel vs. XP Technology**

We would like to cut through the competitors claims as well as their edges.



**1-800-551-5172  
to order today!**



American Eagle Instruments® Inc.  
6575 Butler Creek Rd. Missoula, MT 59808  
1-800-551-5172 • FAX 1-406-549-7452 • www.am-eagle.com

**FREEDOM TO CHOOSE!**

# Quik-Tip™

Cone Sockets With XP Technology™

The Only Sharpen-free  
Cone Sockets On the Planet!






American Eagle Instruments® Inc.

# FREEDOM TO CHOOSE NOT TO SHARPEN!

## Quik-Tip™ Cone Sockets with XP Technology!

Our new line of XP Technology™ instruments now include the convenience of cone socket re-tipping. We continue our commitment to bring our customers the best instruments. For more information on XP Technology please visit us at [www.am-eagle.com/xp](http://www.am-eagle.com/xp) or call 1-800-551-5172.



## HANDLE OPTIONS



AEQTHX - 3/8" - 9.5mm Resin Cone Socket Handle



AEQTHZ - 3/8" - 9.5mm S.S. Cone Socket Handle

## ACCESSORIES



AEQTW - 1/8" - S.S. Cone Socket Wrench